TO DETERMINE YOUNG'S MODULUS OF ELASTICITY OF THE MATERIAL OF A BAR BY THE METHOD OF FLEXURE

Aim: To determine the Young's modulus of elasticity of the given material.

Apparatus: meter scale, weights, screw gauge, spherometer and slide caliper

Theory: If a low-weight bar of breadth b and depth d is placed horizontally on two knife-edges separated by a distance L, and a load of mass m, applied at the mid-point of the bar, produces a depression l of the bar, then Young's modulus Y of the material of the bar is given by

$$Y = \frac{gL^3}{4bd^3} \frac{m}{l} \tag{1}$$

where g is the acceleration due to gravity. This is the working formula of the experiment, and is valid so long as the slope of the bar at any point with respect to the unstrained position is much less than unity. Here Y is determined by measuring the quantities b, d, L and the mean depression l corresponding to a load m. If b, d, L and l are measured in cm, m in gm, g is expressed in cm/sec², and then Y is obtained in dyne/cm².

Procedure:

- (i) Measure the length of the given bar with a meter scale and place the centre of the bar at the middle point of the two supports fixed to the table.
- (ii) Place rectangular hook cum hanger at the centre of the scale. Place the spherometer such that the screw of the spherometer is on top of the rectangular hook cum hanger.
- (iii) Make the electrical connections as shown in the Figure (1) by connecting the power supply, spherometer and galvanometer.
- (iv) Switch on the power supply, adjust the circular scale such that the screw of the spherometer (S) touches the rectangular hook, indicating deflection in the Galvanometer (G).
- (v) Take the initial reading on the spherometer.
- (vi) Hang a weight to the hook and then you notice that galvanometer reading come back to zero (why?). Again adjust the spherometer such that screw touches the hook and galvanometer deflects. Take the depression of the bar reading corresponding the weight hung.
- (vii) Continue the same procedure for all the weights. Avoid back lash error while using the spherometer.
- (viii) Repeat the experiment once again with all the weights.



Figure 1: Schematic diagram of electrical connections

- (ix) Determine the vernier constant (least count) of the slide callipers and measure with it the breadth b of the bar at three different places. Calculate the mean breadth of the bar. Note the zero error, if any, of the slide callipers and find the correct value of b.
- (x) Determine the least count of the screw gauge and measure the depth d of the bar at a number of places along the length of the bar. Find the mean value. Note the zero error, if any, of the screw gauge and obtain the correct value of d.
- (xi) Draw a graph with the load m in grams along the X-axis and the corresponding depression l in cm along the Y-axis and determine the value of Y.

Experimental Results:

Table-1

Least count of spherometer =..... cm

Pitch of the screw p (cm)	No. of divisions n on the	Least count = p/n (cm)
	circular scale	

Table-2 :Load-depression data for chosen length

No.	Load	spherometer reading for		Spherometer reading for			Mean	Depress	
of	in	increasing load (cm)-			increas	sing load	(cm)-	reading	ion
obs.	(kg)	first measurement			second measurement			(cm)	l(cm)
		Main	Circular	Total	Main	Circular	Total		
		scale			scale				
1								(a)	0
2								(b)	(b) - (a)
3								(c)	(c) - (a)
••		••	••			••	••	•••	•••

Distance between the knife-edges $L = \dots \dots cm$

To get the average $\frac{m}{l}$, plot load (m) in kg versus Depression (l) and obtain slope 's' and error in slope ' Δ s' by straight line fit.

Table-3

Least count (l.c.) of the slide caliper

..... Divisions of the Vernier scale = divisions of the main scale.

Value of smallest main	Value of smallest Vernier scale	Least count
scale division (l_1)	division $(l_2 = \frac{n}{m}l_1)$ (cm)	$(l.c.) = (l_1 - l_2) (cm)$

Table-4: Measurement of breadth (*b*) of the bar by slide caliper

No.	Readings	(cm) of	Total	Mean b	Zero error	Corrected b
of	the		reading	(cm)	(cm)	(cm)
obs.	Main	Vernier	<i>b</i> (cm)			
	scale					

Table-5: Least count (l.c.) of the screw gauge

No. of divisions n on the	Least count = p/n (cm)
circular scale	
	•••••
	No. of divisions n on the circular scale

Table–6: Measurement of depth (*d*) of the bar by the screw gauge

No.	Readings (cm)	of the	Total	Mean d	Zero error	Correct d
of	Main scale Vernier		reading	(cm)	(cm)	(cm)
obs.			d (cm)			
1	•••					
2						
3						

Discussion:

Even though the entire scale or bar deforms by applying the stress to given bar or scale, length L in $Y = \frac{gL^3}{4bd^3} \frac{m}{l}$ (or $Y = \frac{gL^3}{4bd^3} s$) corresponds to the portion between the two knife edges but not the total length of the given bar or meter scale. *Why?* Refer to supporting material.

Computation of error: The error in Y due to propagation of errors in the measurement of L, b, d and s is given by

$$\frac{\Delta Y}{Y} = \sqrt{\left(3\frac{\Delta L}{L}\right)^2 + \left(\frac{\Delta b}{b}\right)^2 + \left(3\frac{\Delta d}{d}\right)^2 + \left(\frac{\Delta s}{s}\right)^2}$$

Precautions: (i) In the expression for Y, both the length L between the knife-edges and the depth d of the bar occur in powers of three. But as d is much smaller than L, much care should be taken while measuring d in order to minimize the error in Y.

- (ii) Care should be taken to make the beam horizontal and to load the bar at its mid-point.
- (iii) Try to avoid parallax and back-lash errors during measurements.

Questions: How does the weight of the bar affect the experiment?